

TRUE MANAGER AIR

Truth Behind The Numbers

Track your data to get valuable insights that link daily activity with results. Share data with healthcare providers, friends and family to help manage diabetes.

The new TRUE MANAGER™ AIR App features:

- Simple to navigate logbook that shows date, time, blood glucose results and event tags from meter
- Each blood glucose result is color coded to indicate if within acceptable ADA recommended ranges* for people with diabetes

TRUE MANAGER™ AIR – Insight At Your Fingertips



Upgrade to the TRUE METRIX® AIR Self-Monitoring Blood Glucose System



TRUE
METRIX
AIR

TECHNOLOGY WITH REACH

PERFORMANCE | KNOWLEDGE | TECHNOLOGY | ACCURACY



TRUE
METRIX
AIR

†Go to www.trividiahealth.com/truemanagerair for a list of compatible mobile devices and applications. ‡ American Diabetes Association. Diabetes advocacy. Sec. 14. In Standards of Medical Care in Diabetes 2015. Diabetes Care 2015;38(Suppl. 1):S86–S87; Consult with your doctor or healthcare professional for your target ranges and how often you should test.

NICO-3622HMT 01/18 © 2018 TRUE METRIX, TRIPLE SENSE TECHNOLOGY and TRUE MANAGER, are trademarks of Trividia Health, Inc. Health Mart, the Health Mart logo and the Health Mart Private label logo are registered trademarks of Health Mart Systems, Inc. All rights reserved.

PERFORMANCE IS CONFIDENCE

TRUE METRIX® AIR Advanced Glucose Management

With integrated wireless capability and the TRUE MANAGER™ AIR App, results are seamlessly delivered to your mobile device† to provide simple tracking and insights.



FEATURES:

- No Coding
- Tiny 0.5 microliter sample size
- As fast as 4 seconds
- Stores 1,000 results
- Data management
- Download capabilities

TECHNOLOGY IS CONFIDENCE



TRUE METRIX® Test Strips – Featuring TRIPLE SENSE TECHNOLOGY®

Advanced technologies – the meter, a complex algorithm, chemistry and electrodes - on the test strip, work together as part of the TRUE METRIX® AIR system to produce accurate results.

TRIPLE SENSE TECHNOLOGY® provides proven accuracy and confidence in results.



The TRUE METRIX® AIR Test Strips are designed to detect, analyze and correct for common variables such as hematocrit, which is the number of red blood cells in your blood sample, and temperature, which could impact the accuracy of results.

KNOWLEDGE IS CONFIDENCE

EVENT TAGGING

Easy-to-use TRUE METRIX® AIR with advanced event tags empowers you to make the connection between personal lifestyle and results.



BEFORE MEAL

AFTER MEAL

EXERCISE

SICK

MEDICATION

OTHER


ACCURACY IS CONFIDENCE

Clinical study results concluded that first-time users are able to achieve accurate results when using the TRUE METRIX® AIR system. Additionally, users rated the TRUE METRIX® AIR system as easy-to-use and rated the instructions for use as clear and easy to understand.

Together, these results indicate that the TRUE METRIX® AIR system is accurate and easy-to-use for people with diabetes.¹

Customers Give High Ratings to the TRUE METRIX® AIR System

Evaluation of Ease-of-Use With TRUE METRIX® AIR



Responses on a scale of 1 to 5, with 1 = strongly disagree and 5 = total agreement.

1. Data on file.